
```
# IMPORT PACKAGES CELL
# python scientific computing
import numpy as np
# python plotting library
import matplotlib.pyplot as plt
import matplotlib.colors as colors
#from matplotlib.patches import Rectangle
%matplotlib inline
# python curve fitting library
from scipy.optimize import curve fit
from scipy.stats import norm
# matplotlib styles for scientific figures
# https://github.com/garrettj403/SciencePlots
# https://pypi.org/project/SciencePlots/
import scienceplots
plt.style.use(['science', 'grid']) # 'science', 'ieee', 'grid'
# default figure size
plt.rcParams['figure.figsize'] = [8, 6]
# general font size (affects all text)
plt.rcParams['font.size'] = 16
# axis label font size
plt.rcParams['axes.labelsize'] = 16
# x-axis tick label size
plt.rcParams['xtick.labelsize'] = 14
# y-axis tick label size
plt.rcParams['ytick.labelsize'] = 14
# title font size
plt.rcParams['axes.titlesize'] = 18
# change size of plots
plt.rcParams['figure.dpi'] = 100
# SAVE PLOT
# dpi resolution to save images
image save dpi = 1200
path='graphs/'
format='.pdf'
save = 0
```

Depletion voltage

```
# Open the file in read mode
file = open('Current-Voltage-Char.txt', 'r')
line = file.readline()
line = file.readline()
```

```
voltage = []
current = []
while line:
    #print(line.split())
    voltage.append(float(line.split()[0])) # Extract voltage
    current.append(float(line.split()[2])) # Extract current
    #print(voltage[-1], current[-1])
    # Read the next line
    line = file.readline() # Read the next line
# Close the file
file.close()
voltage = np.array(voltage)
current = np.array(current)
plt.figure()
plt.plot(voltage, current,
         marker='o',
         markerfacecolor=(0,0,1,0.5),
         markeredgecolor=(0,0,1,1.0),
         markeredgewidth=1.2,
         lw=1.
         ls='none',
         markersize=12,
         label='Leakage current')
plt.axvline(60, color='r', lw=2, ls='--',
            label='Depletion voltage')
plt.title('Current-voltage characteristic', pad=20)
plt.xlabel('Voltage [V]')
plt.ylabel('Current [$\\mu$A]')
plt.xlim(0,)
plt.ylim(0.7, 2)
plt.xticks(list(plt.xticks()[0]) + [60])
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='lower right')
plt.tight layout()
if save==1: plt.savefig(path+'CurrVoltChar'+format,
                        dpi=image_save_dpi, bbox inches='tight')
plt.show()
```

Current-voltage characteristic



Pedestal Noise

```
file = open('Pedestal.txt', 'r')
adc = np.loadtxt(file, delimiter=';')
file.close()
pedestal = np.zeros(128)
shift = np.zeros(1000)
noise = np.zeros(128)


for i in range(128):
    pedestal[i] = np.mean(adc[i, :])

for i in range(1000):
    shift[i] = np.mean(adc[:, i] - pedestal[:])

for i in range(128):
    noise[i] = np.sqrt( (1/999)*np.sum((adc[i, :] - pedestal[i] - shift[:])**2) )
```


```
xbins = np.arange(1, 129)
plt.figure()
plt.bar(xbins, pedestal,
        width=1,
        facecolor=(0,0,1,0.5),
        edgecolor=(0,0,1,1.0),
        lw=1,
        ls='solid') #,
        #label='LABEL')
plt.title('Pedestal',pad=20)
plt.xlabel('Strip number')
plt.ylabel('Entries')
plt.xlim(-1, 130)
plt.ylim(507, 518)
plt.ticklabel format(style='sci', axis='both', scilimits=[-4,4])
#plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'pedestal'+format,
                        dpi=image_save_dpi, bbox_inches='tight')
plt.show()
```

Pedestal




```
plt.figure()
plt.bar(xbins, noise,
        width=1,
        facecolor=(0,0,1,0.5),
        edgecolor=(0,0,1,1.0),
        lw=1,
        ls='solid') #,
        #label='LABEL')
plt.title('Noise',pad=20)
plt.xlabel('Strip number')
plt.ylabel('Entries')
plt.xlim(-1, 130)
plt.ylim(1.6, 2.3)
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
\#plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'noise'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Noise



```
hist, bin edges = np.histogram(shift, bins=20, density=True)
def gaussian(x, amp, mu, sigma):
    """Gaussian function for fitting"""
    return amp * np.exp(-(x - mu)**2 / (2 * sigma**2))
bin centers = (bin edges[:-1] + bin edges[1:]) / 2
xgauss = np.linspace(-10, 10, 1000)
# Fit Gaussian to the distribution
params, cov = curve fit(gaussian, bin_centers, hist)
amp, mu, sigma = params
fit_curve = gaussian(xgauss, amp, mu, sigma)
plt.figure()
plt.hist(shift,
         bins=20,
         density=True,
         histtype='barstacked',
         facecolor=(0,0,1,0.5),
         edgecolor=(0,0,1,1.0),
         lw=1,
         ls='solid',
         label='Data')
plt.plot(xgauss, fit_curve, 'r-', linewidth=2,
         label=f'Gaussian Fit:\n $\\mu$ = {mu:.2f}, $\\sigma$ =
{sigma:.2f}')
plt.title('Common mode shift / ADC counts', pad=20)
plt.xlabel('Shift')
plt.ylabel('Entries (normalized)')
plt.xlim(-9, 9)
plt.ylim(0, )
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'shift'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Common mode shift / ADC counts



Calibration measurements

Delay measurement

```
file = open('Calib/delay_optimum_measure.txt', 'r')

next(file)
next(file)
line = file.readline()

time = []
adc_counts = []

while line:
    time.append(line.split()[0])
    adc_counts.append(line.split()[1])


line = file.readline()

file.close()

# Replace commas with periods and convert to float
time_float = np.array([float(t.replace(',', '.'))) for t in time])
```

```
adc counts float = np.array([float(t.replace(',', '.')) for t in
adc counts])
plt.figure()
plt.plot(time float, adc counts float,
         marker='o',
         markerfacecolor=(0,0,1,0.5),
         markeredgecolor=(0,0,1,1.0),
         markeredgewidth=1.2,
         lw=1,
         ls='none',
         markersize=12,
         label='Measured counts')
plt.axvline(64, color='r', lw=2, ls='--',
            label='Optimal delay value')
plt.title('Delay measurement', pad=20)
plt.xlabel('Time [ns]')
plt.ylabel('ADC counts')
plt.xlim(0, 250)
plt.ylim(0, )
plt.xticks(list(plt.xticks()[0]) + [64])
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'optDelay'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Delay measurement



Channel measurement

```
def data_read(name):
    Reads the calibration data from a file and returns the calibration
values
    """

    file = open(name+'.txt', 'r')
    next(file)
    next(file)
    line = file.readline()


    charge = []
    adc = []

while line:
        charge.append(line.split()[0])
        adc.append(line.split()[1])
        line = file.readline()


file.close()
```

```
# Replace commas with periods and convert to float
    charge float = np.array([float(t.replace(',', '.')) for t in
charge])
    adc float = np.array([float(t.replace(',', '.')) for t in adc])
    return charge_float, adc_float
def plot channel(cha, adc, color, linest, cha name):
    plt.plot(cha, adc,
              marker='o',
              markerfacecolor=(color, 0.5),
              markeredgecolor=(color, 0.8),
              markeredgewidth=1,
              lw=2.
              ls=linest,
              color=color,
              markersize=6,
              label=cha name)
    return
cha20, adc20 = data read('Calib/ch20 plot')
cha40, adc40 = data read('Calib/ch40 plot')
cha60, adc60 = data read('Calib/ch60 plot')
cha80, adc80 = data read('Calib/ch80 plot')
cha80_0V, adc80_0V = data_read('Calib/ch80_0_VOLT_plot')
cha100, adc100 = data read('Calib/ch100 plot')
average = (adc20 + adc40 + adc60 + adc80 + adc100) / 5
color = [[1,0,0], [0,1,0], [0,0,1], [1,0,1], [0,1,1]]
plt.figure()
plot_channel(cha20, adc20, color[0], 'none', 'Channel 20')
plot_channel(cha40, adc40, color[1], 'none', 'Channel 40')
plot_channel(cha60, adc60, color[2], 'none', 'Channel 60')
plot_channel(cha80, adc80, color[3], 'none', 'Channel 80')
plot channel(cha100, adc100, color[4], 'none', 'Channel 100')
plt.plot(cha20, average,
          color='black',
          lw=2,
          ls='-'
          label='Mean value')
plt.title('Calibration curve (5 channels)', pad=20)
plt.xlabel('Injected charge [e]')
plt.ylabel('ADC counts')
```

Calibration curve (5 channels)



Calibration curve (0 Volt)



Polynomial fit


```
def polyfun(x, a, b, c, d, e):
    """4-th degree polynomial fitting"""
    return a*x**4 + b*x**3 + c*x**2 + d*x + e

xspace = np.linspace(0, 300, 1000)

params, cov = curve_fit(polyfun, average[0:120], cha20[0:120])
a, b, c, d, e = params
fit_curve = polyfun(xspace, a, b, c, d, e)
```

```
plt.figure()
plot channel(average, cha20, color[0], 'none', 'Mean ADC counts')
plt.plot(xspace, fit curve,
         color='blue',
         lw=2,
         ls='-',
         label='Fit')
plt.axvline(250, color='k', lw=2, ls='--',
            label='Cut of the fit range')
plt.title('Fit of 4-th degree polynomial to the mean'
          '\n'
          'ADC counts of the calibration run', pad=20)
plt.xlabel('Mean ADC counts')
plt.ylabel('Injected charge [e]')
plt.xlim(0, 280)
plt.ylim(0, )
plt.ticklabel format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'polyFit'+format,
                        dpi=image_save_dpi, bbox_inches='tight')
plt.show()
```

Fit of 4-th degree polynomial to the mean ADC counts of the calibration run




```
# Print the parameters of the polynomial fit
# using scientific notation {a:e} as example
print(f'a = {a:e}')
print(f'b = {b:e}')
print(f'c = {c:e}')
print(f'd = {d:e}')
print(f'e = {e:e}')

a = 5.263758e-05
b = -1.913336e-02
c = 2.693738e+00
d = 1.753607e+02
e = 9.804066e+02
```

Measuring the strip sensor by using the laser


```
delay, adc_laser = data_read('Laserscan/laser_opt_plot')
plt.figure()
plot_channel(delay, adc_laser, color[0], 'none', 'ADC counts')
```

Results of the Laser Sync run



```
lascan = np.loadtxt('Laserscan.txt',unpack=True)
plt.figure()
plt.imshow(lascan, vmin=0, vmax=100, cmap='viridis', origin='lower',
```


Results of the Laser Sync run



```
plt.figure()


plt.imshow(lascan, vmin=0, vmax=100, cmap='viridis', origin='lower',
aspect='auto')
plt.colorbar(label='ADC counts')
```

Results of the Laser Sync run




```
label='Start of peak')
plt.axvline(3, color='g', lw=2, ls='-.',
             label='Maximum')
plt.axvline(7, color='r', lw=2, ls='--',
             label='End of peak')
plt.axvline(12, color='k', lw=2, ls='--')
plt.axvline(15, color='g', lw=2, ls='-.')
plt.axvline(19, color='r', lw=2, ls='--')
plt.title('Strip\'s width and laser extension (Channel 73)', pad=20)
plt.xlabel('Distance / 10 [$\\mu$m]')
plt.ylabel('ADC counts')
plt.xlim(-5, 35)
plt.ylim(0, )
plt.xticks(list([0, 3, 7, 12, 15, 19, 25, 30, 35]))
plt.ticklabel format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight_layout()
if save==1: plt.savefig(path+'width73'+format,
                          dpi=image_save_dpi, bbox_inches='tight')
plt.show()
```

Strip's width and laser extension (Channel 73)



```
plt.figure()
plot_channel(np.linspace(0, 34, 35), lascan[73, :],
              color[0], '--', 'Channel 74')
plt.axvline(15, color='k', lw=2, ls='--',
             label='Start of peak')
plt.axvline(19, color='g', lw=2, ls='-.',
             label='Maximum')
plt.axvline(23, color='b', lw=2, ls='--',
            label='End of peak')
plt.axvline(28, color='k', lw=2, ls='--')
plt.axvline(31, color='g', lw=2, ls='-.')
plt.axvline(34, color='b', lw=2, ls='--')
plt.title('Strip\'s width and laser extension (Channel 74)', pad=20)
plt.xlabel('Distance / 10 [$\\mu$m]')
plt.ylabel('ADC counts')
plt.xlim(0, 35)
plt.ylim(0, )
plt.xticks(list([5, 10, 15, 19, 23, 28, 31, 34]))
plt.ticklabel format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'width74'+format,
                         dpi=image save dpi, bbox inches='tight')
plt.show()
```

Strip's width and laser extension (Channel 74)



```
1.1.1
import os
def combine text files(input files, output file, delimiter='\t'):
    Combines multiple single-column text files into a multi-column
text file.
    Args:
        input files (list): List of paths to input text files
        output_file (str): Path for the output combined text file
        delimiter (str): Column separator for output file (default:
tab)
    # Read all input files
    columns = [1]
    for file_path in input files:
        with open(file_path, 'r') as f:
            # Read all lines and strip whitespace
            column = [line.strip() for line in f.readlines()]
            columns.append(column)
```


```
# Verify consistent row count
    row counts = [len(col) for col in columns]
    if len(set(row counts)) > 1:
        raise ValueError(f"Inconsistent row counts: {row counts}. All
files must have same number of rows.")
    # Transpose columns into rows
    combined data = []
    for i in range(len(columns[0])):
        row = delimiter.join(col[i] for col in columns)
        combined data.append(row)
    # Write output file
    with open(output file, 'w') as f:
        f.write('\n'.join(combined data))
'\nimport os\n\ndef combine text files(input files, output file,
                        """√n
delimiter=\'\t\'):\n
                                 Combines multiple single-column text
files into a multi-column text file.\n\n
                                                            input files
                                            Args:\n
(list): List of paths to input text files\n
                                                   output file (str):
Path for the output combined text file\n
                                                delimiter (str):
                                                     """\n
Column separator for output file (default: tab)\n
                                                              # Read
all input files\n
                     columns = []\n for file path in input files:\n
with open(file path, \'r\') as f:\n
                                               # Read all lines and
                              column = [line.strip() for line in
strip whitespace\n
f.readlines()]\n
                            columns.append(column)\n\n
                                                          # Verify
consistent row count\n
                          row counts = [len(col) for col in columns]\n
if len(set(row counts)) > 1:\n
                                      raise ValueError(f"Inconsistent
row counts: {row counts}. All files must have same number of rows.")\
       # Transpose columns into rows\n
                                          combined data = []\n
n\n
                                      row = delimiter.join(col[i] for
i in range(len(columns[0])):\n
                         combined data.append(row)\n\n
col in columns)\n
                with open(output file, \'w\') as f:\n
output file\n
f.write(\'\n\'.join(combined data))\n'
combine text files(
    input files=['CCEL/0CCEL.txt',
                 'CCEL/10CCEL.txt',
                 'CCEL/20CCEL.txt'
                 'CCEL/30CCEL.txt'
                 'CCEL/40CCEL.txt'
                 'CCEL/50CCEL.txt'
                 'CCEL/60CCEL.txt',
                 'CCEL/70CCEL.txt'
                 'CCEL/80CCEL.txt',
                 'CCEL/90CCEL.txt'
                 'CCEL/100CCEL.txt'
                 'CCEL/110CCEL.txt',
```

```
'CCEL/120CCEL.txt',
                 'CCEL/130CCEL.txt'
                 'CCEL/140CCEL.txt'
                 'CCEL/150CCEL.txt'
                 'CCEL/160CCEL.txt'
                 'CCEL/170CCEL.txt'
                 'CCEL/180CCEL.txt'
                 'CCEL/190CCEL.txt'
                 'CCEL/200CCEL.txt',
   output file='CCEL combined.txt',
   delimiter='\t'
)
                            input files=['CCEL/0CCEL.txt',\n
"\ncombine text files(\n
'CCEL/10CCEL.txt',\n
                                      'CCEL/20CCEL.txt',\n
'CCEL/30CCEL.txt',\n
                                      'CCEL/40CCEL.txt',\n
'CCEL/50CCEL.txt',\n
                                      'CCEL/60CCEL.txt',\n
'CCEL/70CCEL.txt',\n
                                      'CCEL/80CCEL.txt',\n
                                      'CCEL/100CCEL.txt',\n
'CCEL/90CCEL.txt',\n
'CCEL/110CCEL.txt',\n
                                       'CCEL/120CCEL.txt',\n
'CCEL/130CCEL.txt',\n
                                       'CCEL/140CCEL.txt',\n
'CCEL/150CCEL.txt',\n
                                       'CCEL/160CCEL.txt',\n
'CCEL/170CCEL.txt',\n
                                       'CCEL/180CCEL.txt',\n
'CCEL/190CCEL.txt',\n
                                       'CCEL/200CCEL.txt',\n
        output file='CCEL combined.txt',\n
                                               delimiter='\t'\n)\n"
],\n
```


Charge collection efficiency

CCE Laser

Heatmap of the ADC counts per each channel and bias voltage



Heatmap of the ADC counts per each channel and bias voltage



```
def deepFun(x, a):
    D = 300
    U_dep = 80
    return (1 - np.exp(-(D * np.sqrt(x/U_dep))/a)) / (1 - np.exp(-D/a))


xspace = np.linspace(0, 80, 1000)

data = ccel[71, :]/ccel[71, 8]

params, cov = curve_fit(deepFun, np.linspace(0, 80, 9), data[:9],
```

```
bounds=(0, 300))
a opt = params[0]
fit curve = deepFun(xspace, a opt)
a err = np.sqrt(cov[0, 0])
print(f"Optimal a parameter: {a opt:.6f}")
print(f"Uncertainty in a: ± {a_err:.6f}")
print(f"Relative uncertainty: ± {100*a err/a opt:.4f}%")
Optimal a parameter: 252.717907
Uncertainty in a: ± 30.108226
Relative uncertainty: ± 11.9138%
plt.figure()
plot channel(np.linspace(0, 200, 21), data,
             color[0], '--', 'Channel 72')
plt.axvline(80, color='k', lw=2, ls='--',
            label='Start of plateau $U {dep}$')
plt.plot(xspace, fit curve,
         color='blue',
         lw=2,
         ls='-'
         label='Fit')
plt.title('Charge collection efficiency (Channel 72)', pad=20)
plt.xlabel('Voltage [V]')
plt.ylabel('CCE (normalized)')
plt.xlim(-5, 205)
plt.ylim(0, )
plt.xticks(np.arange(0, 201, 10))
plt.xticks(rotation=90)
plt.xticks(list(plt.xticks()[0]) + [80])
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight lavout()
if save==1: plt.savefig(path+'channel72'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Charge collection efficiency (Channel 72)



CCE Beta source

```
import os

def source_avg(input_files, output_file):
    Compute first column averages for multiple files and save to single output file

Args:
    input_files (list): List of input file paths output_file (str): Path for output text file

"""

averages = []

for file_path in input_files:
    values = []


with open(file_path, 'r') as f:
    for line in f:
    # Skip empty lines and lines starting with #
```

```
if not line.strip() or line.startswith('#'):
                    continue
                # Split line into columns
                columns = line.split()
                # Process only if there's at least one column
                if columns:
                    try:
                        # Extract first column value and convert to
float
                        value = float(columns[0])
                        values.append(value)
                    except ValueError:
                        # Skip lines with non-numeric values in first
column
                        continue
        # Calculate average if we have valid values
        if values:
            avg = sum(values) / len(values)
            averages.append(avg)
        else:
            # Use 0 if no valid values found
            averages.append(0.0)
    # Write all averages to output file
    with open(output file, 'w') as f:
        for avg in averages:
            f.write(f"{avg}\n")
   print(f"Saved {len(averages)} averages to {output file}")
1.1.1
'\nimport os\n\ndef source avg(input files, output file):\n """\n
Compute first column averages for multiple files and save to single
output file\n\n
                   Args:\n
                                  input files (list): List of input
                    output file (str): Path for output text file\n
file paths\n
         averages = []\n\n for file path in input files:\n
                       with open(file_path, \'r\') as f:\n
values = []\n\n
for line in f:\n
                                # Skip empty lines and lines starting
with #\n
                        if not line.strip() or
line.startswith(\'#\'):\n
                                             continue\n\n
                                           columns = line.split()\n\n
# Split line into columns\n
# Process only if there\'s at least one column\n
                                                                 if
columns:\n
                                                            # Extract
                              trv:\n
first column value and convert to float\n
                                                                  value
= float(columns[0])\n
                                             values.append(value)\n
except ValueError:\n
                                            # Skip lines with non-
numeric values in first column\n
                                                        continue\n\n
```

```
# Calculate average if we have valid values\n
                                                      if values:\n
avg = sum(values) / len(values)\n
                                              averages.append(avg)\n
else:\n
                   # Use 0 if no valid values found\n
averages.append(0.0)\n\n
                            # Write all averages to output file\n
with open(output file, \'w\') as f:\n
                                              for avg in averages:\n
f.write(f''\{avg\}\setminus n'')\setminus n\setminus print(f''Saved \{len(averages)\} averages to
{output file}")\n'
source avg(
    input files=['source measurement/0 Cluster adc entries.txt',
                  'source measurement/10 Cluster adc entries.txt'
                 'source measurement/20 Cluster adc entries.txt'
                  'source measurement/30 Cluster adc entries.txt'
                 'source measurement/40 Cluster adc entries.txt'
                 'source measurement/50 Cluster adc entries.txt'
                 'source measurement/60 Cluster adc entries.txt'
                 'source measurement/70 Cluster adc entries.txt'
                 'source measurement/80 Cluster adc entries.txt'
                 'source measurement/90 Cluster adc entries.txt'
                 'source measurement/100 Cluster adc entries.txt'
                 'source measurement/110 Cluster adc entries.txt'
                 'source measurement/120 Cluster adc entries.txt'
                 'source measurement/130 Cluster adc entries.txt'
                 'source measurement/140 Cluster adc entries.txt'
                 'source_measurement/150_Cluster_adc_entries.txt'
                 'source measurement/160 Cluster adc entries.txt'
                 'source measurement/170 Cluster adc entries.txt'
                 'source measurement/180 Cluster adc entries.txt'
                 'source measurement/190 Cluster adc entries.txt'
                 'source measurement/200 Cluster adc entries.txt',
    output file='SOURCE combined.txt'
)
"\nsource avg(\n
input files=['source measurement/0 Cluster adc entries.txt',\n
'source measurement/10 Cluster adc_entries.txt',\n
'source measurement/20 Cluster adc entries.txt',\n
'source measurement/30 Cluster adc entries.txt',\n
'source measurement/40 Cluster adc entries.txt',\n
'source measurement/50 Cluster adc entries.txt',\n
'source_measurement/60_Cluster_adc entries.txt',\n
'source measurement/70 Cluster adc entries.txt',\n
'source_measurement/80_Cluster_adc_entries.txt',\n
'source measurement/90 Cluster adc entries.txt',\n
'source measurement/100 Cluster adc entries.txt',\n
'source measurement/110 Cluster adc entries.txt',\n
'source measurement/120 Cluster adc entries.txt',\n
```


```
'source measurement/130 Cluster adc entries.txt',\n
'source measurement/140 Cluster adc entries.txt',\n
'source measurement/150 Cluster adc entries.txt',\n
'source measurement/160 Cluster adc entries.txt',\n
'source measurement/170 Cluster adc entries.txt',\n
'source_measurement/180_Cluster_adc_entries.txt',\n
'source measurement/190 Cluster adc entries.txt',\n
'source measurement/200 Cluster adc entries.txt',\
                  ],\n output file='SOURCE combined.txt'\n)\n"
file = open('SOURCE combined.txt', 'r')
line = file.readline()
src = []
while line:
    src.append(line.split()[0])
    line = file.readline()
file.close()
src float = np.array([float(t.replace(',', '.')) for t in src])
plt.figure()
plot_channel(np.linspace(0, 200, 21), src_float/src_float[8],
             color[2], '--', 'Beta source')
plt.axvline(80, color='k', lw=2, ls='--',
            label='Start of plateau $U {dep}$')
plt.title('Charge collection efficiency (Beta source)', pad=20)
plt.xlabel('Voltage [V]')
plt.ylabel('CCE (normalized)')
plt.xlim(-5, 205)
plt.ylim(0, )
plt.xticks(np.arange(0, 201, 10))
plt.xticks(rotation=90)
plt.xticks(list(plt.xticks()[0]) + [80])
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'betaSource'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Charge collection efficiency (Beta source)




```
plt.figure()
plot channel(np.linspace(0, 200, 21), data,
             color[0], '--', 'Channel 72 (Laser)')
plot_channel(np.linspace(0, 200, 21), src_float/src_float[8],
             color[2], '--', 'Beta source')
plt.axvline(80, color='k', lw=2, ls='--',
            label='Start of plateau $U_{dep}$')
plt.title('Charge collection efficiency (Beta source)', pad=20)
plt.xlabel('Voltage [V]')
plt.ylabel('CCE (normalized)')
plt.xlim(-5, 205)
plt.ylim(0, )
plt.xticks(np.arange(0, 201, 10))
plt.xticks(rotation=90)
plt.xticks(list(plt.xticks()[0]) + [80])
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
plt.legend(loc='best')
plt.tight layout()
```

Charge collection efficiency (Beta source)




Large source scan

```
(bins[i], cluster[i]), # Anchor point (top of
bar)
                     xytext=(0, 3),
                                                 # 3-point vertical
offset
                     textcoords='offset points', # Offset in points
                     ha='center',
                                                 # Horizontal alignment
                     va='bottom',
                                                 # Vertical alignment
                      rotation=0,
                      fontsize=12)
#plt.title('Number of clusters', pad=20)
plt.xlabel('Number of clusters')
plt.ylabel('Entries')
plt.xlim(-0.5,)
#plt.ylim(0, )
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
#plt.legend(loc='best')
plt.tight_layout()
if save==1: plt.savefig(path+'cluster'+format,
                         dpi=image save dpi, bbox inches='tight')
plt.show()
```




```
cluchan = np.array([0,
                    562927.
                    406888,
                    47034.
                    10298,
                    3375,
                    1364,
                    566,
                    262,
                    137,
                    77,
                    30,
                    12,
                    7,
                    5,
                    3,
                    1,
                    0,
                    0])
clubins = np.arange(len(cluchan))
plt.figure()
plt.bar(clubins, cluchan, #color='blue', edgecolor='black')
        width=1.
        facecolor=(0,0,1,0.5),
        edgecolor=(0,0,1,1.0),
        lw=1.
        ls='solid') #,
        #label='LABEL')
for i in range(len(cluchan)):
    if cluchan[i] > 0:
        plt.annotate(f'{cluchan[i]:.0f}',
                     (clubins[i], cluchan[i]), # Anchor point (top of
bar)
                     xytext=(0, 3),
                                                # 3-point vertical
offset
                     textcoords='offset points', # Offset in points
                     ha='center',
                                               # Horizontal alignment
                     va='bottom',
                                               # Vertical alignment
                     rotation=90,
                     fontsize=12)
#plt.title('Number of clusters', pad=20)
plt.xlabel('Number of channels per clusters')
plt.ylabel('Entries')
plt.xlim(-0.5,)
plt.ylim(0, 700000)
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
```



```
file = open('source_measurement/hitmap.txt', 'r')
next(file)
line = file.readline()
hit = []
while line:
    hit.append(line.split()[0])
    line = file.readline()
file.close()
hit_float = np.array([float(t.replace(',', '.'))) for t in hit])
```


```
plt.figure()
plt.bar(xbins, hit float, #color='blue', edgecolor='black')
        width=1,
        facecolor=(0,0,1,0.5),
        edgecolor=(0,0,1,1.0),
        lw=1,
        ls='solid') #,
        #label='LABEL')
plt.title('Hitmap displaying the events of each channel', pad=20)
plt.xlabel('Channels')
plt.ylabel('Hits')
plt.xlim(0, 129)
plt.ylim(₀, )
plt.ticklabel_format(style='sci', axis='both', scilimits=[-4,4])
#plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'hitmap'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Hitmap displaying the events of each channel



```
file = open('source measurement/Cluster adc entries.txt', 'r')
next(file)
line = file.readline()
pulse = []
while line:
    pulse.append(line.split()[0])
    line = file.readline()
file.close()
pulse_float = np.array([float(t.replace(',', '.')) for t in pulse])
pulse energy = pulse float / 3.6 # Convert to keV
plt.figure()
h = plt.hist(pulse float,
             bins=200,
             density=False.
             histtype='barstacked',
             facecolor=(0,0,1,0.5),
             edgecolor=(0,0,1,1.0),
             lw=1,
             ls='solid')#,
             #label='Data')
plt.title('Distribution of the ADC counts', pad=20)
plt.xlabel('ADC counts')
plt.ylabel('Number of events')
plt.xlim(min(h[1])-1, max(h[1])+1)
plt.ylim(0, )
plt.ticklabel format(style='sci', axis='both', scilimits=[-4,4])
#plt.legend(loc='best')
plt.tight layout()
if save==1: plt.savefig(path+'ADCdistrib'+format,
                        dpi=image save dpi, bbox inches='tight')
plt.show()
```

Distribution of the ADC counts



```
e_mean = np.mean(pulse_energy)
plt.figure()
h = plt.hist(pulse energy,
             bins=100,
             density=False,
             histtype='step',
             facecolor=(0,0,1,0.5),
             edgecolor=(0,0,1,1.0),
             lw=1,
             ls='solid',
             label='Mesured energy distibution')
# Find the index of the maximum bin count
mpv idx = np.argmax(h[0])
# Calculate the MPV as the center of the corresponding bin
e_mpv = 0.5 * (h[1][mpv_idx] + h[1][mpv_idx + 1])
plt.axvline(e_mean, color='k', lw=2, ls='--',
            label=r'$E {\mathrm{mean}}'+f'= {e mean:.2f}$ [keV]')
```

Distribution of energy

