

```

import pickle
import matplotlib.pyplot as plt
import numpy as np
from scipy.optimize import curve_fit
import pandas as pd
import pathlib as pl
#from uncertainties.unumpy import nominal_values as noms
#from uncertainties.unumpy import std_devs as stds

```

5.1 Spectrometer Measurement

Taking data with the spectrometer to use as background reference

```

# Import the roomlight reference data as a pandas DataFrame
ref_roomlight_on = pd.read_csv("gc_roomlight.txt", sep="\t")
ref_roomlight_off = pd.read_csv("gc_roomlight_off.txt", sep="\t")

```

```
print(ref_roomlight_on, ref_roomlight_off)
```

	#lambda/nm	DC	C1
0	208.295	638.2	641.0
1	208.670	561.0	559.4
2	209.045	528.2	523.2
3	209.419	517.4	517.2
4	209.794	611.4	613.2
...
2043	1129.520	593.2	588.8
2044	1129.890	583.6	587.2
2045	1130.260	583.6	588.8
2046	1130.640	575.4	592.0
2047	1131.010	587.0	584.6

	#lambda/nm	DC	C1
0	208.295	633.0	633.2
1	208.670	551.8	550.2
2	209.045	517.4	517.8
3	209.419	506.0	509.8
4	209.794	600.2	606.4
...
2043	1129.520	543.2	541.6
2044	1129.890	545.2	546.8
2045	1130.260	531.8	535.4
2046	1130.640	536.0	532.2
2047	1131.010	545.6	543.4

```
[2048 rows x 3 columns]
```

```

# Plot the roomlight data with lights on and off
plt.plot(ref_roomlight_on['#lambda/nm'], ref_roomlight_on['DC'],
label='Roomlight ON, Dark Count')

```

```


plt.plot(ref_roomlight_off['#lambda/nm'], ref_roomlight_off['DC'],
label='Roomlight OFF, Dark Count')
plt.plot(ref_roomlight_on['#lambda/nm'], ref_roomlight_on['C1'],
label='Roomlight ON, LED ON',alpha=0.5)
plt.plot(ref_roomlight_off['#lambda/nm'], ref_roomlight_off['C1'],
label='Roomlight OFF, LED ON')
plt.axhline(y=550, color='black', linestyle='--')
plt.axhline(y=55000, color='black', linestyle='--')
plt.text(-100, 550, f' Counts = 550', va='center', ha='left',
color='black')
plt.text(-100, 55000, f' Counts = 55000', va='center', ha='left',
color='black')
plt.xlabel('Wavelength (nm)')
plt.ylabel('Intensity (Counts)')
plt.title('Roomlight Reference Data')
plt.legend()
plt.grid()
plt.yscale('log')
plt.show()


# LED count - Dark Count for roomlight off
ref_roomlight_corr = pd.DataFrame({
    '#lambda/nm': ref_roomlight_off['#lambda/nm'].astype(float),
    'Count_off': (ref_roomlight_off['C1'] -
ref_roomlight_off['DC']).astype(float),
    'Count_on': (ref_roomlight_on['C1'] -
ref_roomlight_on['DC']).astype(float)
})

# Plot the corrected roomlight off data
plt.plot(ref_roomlight_corr['#lambda/nm'],
ref_roomlight_corr['Count_off'], label='Roomlight OFF, LED ON - Dark
Count Corrected', color='violet')
plt.scatter(ref_roomlight_corr['#lambda/nm'],
ref_roomlight_corr['Count_on'], label='Roomlight ON, LED ON - Dark
Count Corrected', color='gold')
plt.axhline(y=0, color='black', linestyle='--')
plt.xlabel('Wavelength (nm)')
plt.ylabel('Intensity (Counts)')
plt.title('Dark Count (Background) Corrected LED Readings')
plt.legend()
plt.grid()
plt.show()

```

Roomlight Reference Data





5.2 Radial Symmetry Measurement

Taking data with varying angles in the spectrometer to test the symmetry of the system

```
# "gc_radial" file path
gc_radial_folder = pl.Path("gc_radial")

# Import Radial Symmetry Dataframe

df_list = []
for file in gc_radial_folder.glob("*.txt"):
    df = pd.read_csv(file, sep="\t", header=None)
    # Add 3rd column: filename without extension
    df['Measurement'] = file.stem
    df.columns = ['#lambda/nm', 'C1', 'Measurement Specifications']
    df_list.append(df[1:])

radial_df = pd.concat(df_list, ignore_index=True)

print(radial_df.head())
```

	#lambda/nm	C1	Measurement Specifications
0	208.295	544.6	DarkCounts_h=10deg_v=34.8deg
1	208.67	487.4	DarkCounts_h=10deg_v=34.8deg

```

2    209.045  453.2  DarkCounts_h=10deg_v=34.8deg
3    209.419  450.4  DarkCounts_h=10deg_v=34.8deg
4    209.794  603.8  DarkCounts_h=10deg_v=34.8deg

# Sum the counts for each measurement specification across all
wavelengths
radial_df['C1'] = radial_df['C1'].astype(float)
summed_radial_df = radial_df.groupby('Measurement Specifications',
as_index=False)[['#lambda/nm', 'C1']].sum(numeric_only=True)

# Drop the '#lambda/nm' column if it STILL exists because it is a
fucking ghost
if '#lambda/nm' in summed_radial_df.columns:
    summed_radial_df = summed_radial_df.drop(columns=['#lambda/nm'])

# Split the first column into three columns
split_cols = summed_radial_df['Measurement
Specifications'].str.split('_', expand=True)

# Extract measurement type, h angle, and v angle
summed_radial_df['Measurement type'] = split_cols[0]
summed_radial_df['h angle'] = split_cols[1].str.extract(r'h=(-\d.]
+deg').astype(float)
summed_radial_df['v angle'] = split_cols[2].str.extract(r'v=(-\d.]
+deg').astype(float)

# Drop the Measurement specifications column
summed_radial_df = summed_radial_df.drop(columns=['Measurement
Specifications'])

# Move the 'C1' column to the end and rename it to 'Counts'
summed_radial_df['Counts'] = summed_radial_df['C1']
summed_radial_df = summed_radial_df.drop(columns=['C1'])

# Subtract DarkCounts from Attenuation for matching h angle and v
angle
Attenuation_df = summed_radial_df[summed_radial_df['Measurement type']
== 'Attenuation']
DarkCount_df = summed_radial_df[summed_radial_df['Measurement type']
== 'DarkCounts']

# Merge Attenuation_df and DarkCount_df on 'h angle' and 'v angle'
subtracted_radial_df = pd.merge(
    Attenuation_df,
    DarkCount_df,
    on=['h angle', 'v angle'],
    suffixes=('_atten', '_dark')
)

# Subtract the counts

```

```
subtracted_radial_df['Counts_subtracted'] =
subtracted_radial_df['Counts_atten'] -
subtracted_radial_df['Counts_dark']
```

```
print(subtracted_radial_df.head())
```

	Measurement	type_atten	h angle	v angle	Counts_atten \
0		Attenuation	-10.0	-2.6	1214415.8
1		Attenuation	-10.0	-6.0	1213050.6
2		Attenuation	-10.0	0.8	1214077.0
3		Attenuation	-10.0	11.0	1208326.8
4		Attenuation	-10.0	14.4	1207325.6

	Measurement	type_dark	Counts_dark	Counts_subtracted
0		DarkCounts	1158286.6	56129.2
1		DarkCounts	1157506.0	55544.6
2		DarkCounts	1158206.0	55871.0
3		DarkCounts	1159245.0	49081.8
4		DarkCounts	1159708.4	47617.2

```
# 2D Histogram of the subtracted radial data
```

```
# Pivot the dataframe to create a 2D grid for plotting
```

```
hist2d_sub = subtracted_radial_df.pivot(index='v angle', columns='h
angle', values='Counts_subtracted')
```

```
plt.figure(figsize=(8, 8))
plt.imshow(hist2d_sub.values,
           extent=[hist2d_sub.columns.min(), hist2d_sub.columns.max(),
                  hist2d_sub.index.min(), hist2d_sub.index.max()],
           origin='lower',
           aspect='auto',
           interpolation='nearest')
plt.colorbar(label='Intensity (Counts)')
```


```
# Find the indices of the maximum value in the histogram
```

```
max_idx = np.unravel_index(np.nanargmax(hist2d_sub.values),
hist2d_sub.shape)
max_v = hist2d_sub.index[max_idx[0]]
max_h = hist2d_sub.columns[max_idx[1]]
plt.axvline(x=max_h, color='red', linestyle='--',
label=f'x={max_h:.1f}')
plt.axhline(y=max_v, color='red', linestyle='--',
label=f'y={max_v:.1f}')
```

```
plt.xlabel('Horizontal angle (deg)')
plt.ylabel('Vertical angle (deg)')
plt.title('2D Histogram of Radial Intensity Measurement at the end of
the Scintillating Fiber with Fixed Excitation Length')
```

```
plt.show()
```

2D Histogram of Radial Intensity Measurement at the end of the Scintillating Fiber with Fixed Excitation Length



5.3 Intensity Measurement

Measure of Intensity with the change of angles (h and v) and fiber excitation distance from the measurement point (x value)

```
# "gc_intensity" file path
gc_intensity_folder = pl.Path("gc_Intensity")

# Import Intensity Dataframe

df2_list = []
```

```

for file in gc_intensity_folder.glob("*.txt"):
    df = pd.read_csv(file, sep="\t", header=None)
    # Add 3rd column: filename without extension
    df['Measurement'] = file.stem
    df.columns = ['#lambda/nm', 'C1', 'Measurement Specifications']
    df2_list.append(df[1:])

```

```

intensity_df = pd.concat(df2_list, ignore_index=True)

```

```

print(intensity_df.head())

```

	#lambda/nm	C1	Measurement Specifications
0	208.295	512.4	Attenuation_h=0deg_v=0deg_x=450mm
1	208.67	469.2	Attenuation_h=0deg_v=0deg_x=450mm
2	209.045	426.4	Attenuation_h=0deg_v=0deg_x=450mm
3	209.419	435	Attenuation_h=0deg_v=0deg_x=450mm
4	209.794	613	Attenuation_h=0deg_v=0deg_x=450mm

```

# Sum the counts for each measurement specification across all
wavelengths

```

```

intensity_df['C1'] = intensity_df['C1'].astype(float)
summed_intensity_df = intensity_df.groupby('Measurement
Specifications', as_index=False)[['#lambda/nm',
'C1']].sum(numeric_only=True)

```

```

# Drop the '#lambda/nm' column if it STILL exists because it is a
fucking ghost

```

```

if '#lambda/nm' in summed_intensity_df.columns:
    summed_intensity_df =
summed_intensity_df.drop(columns=['#lambda/nm'])

```

```

# Split the first column into three columns

```

```

split_cols = summed_intensity_df['Measurement
Specifications'].str.split('_', expand=True)

```

```

# Extract measurement type, h angle, and v angle

```

```

summed_intensity_df['Measurement type'] = split_cols[0]
summed_intensity_df['h angle'] = split_cols[1].str.extract(r'h=([-\.d.]
+deg)').astype(float)
summed_intensity_df['v angle'] = split_cols[2].str.extract(r'v=([-\.d.]
+deg)').astype(float)
summed_intensity_df['x length'] = split_cols[3].str.extract(r'x=([-\.
d.] +)mm').astype(float)

```

```

# Drop the Measurement specifications column

```

```

summed_intensity_df = summed_intensity_df.drop(columns=['Measurement
Specifications'])

```

```

# Move the 'C1' column to the end and rename it to 'Counts'

```

```

summed_intensity_df['Counts'] = summed_intensity_df['C1']

```



```

summed_intensity_df = summed_intensity_df.drop(columns=['C1'])

# Subtract Darkcounts from each Attenuation measurement
Attenuation_intensity_df =
summed_intensity_df[summed_intensity_df['Measurement type'] ==
'Attenuation']
Attenuation_intensity_df =
Attenuation_intensity_df.drop(columns=['Measurement type'])
DarkCount_intensity_df =
summed_intensity_df[summed_intensity_df['Measurement type'] ==
'DarkCounts']
DarkCount_intensity_df =
DarkCount_intensity_df.drop(columns=['Measurement type'])

# Merge Attenuation_df and DarkCount_df on 'h angle', 'v angle' and 'x
length'
subtracted_intensity_df = pd.merge(
    Attenuation_intensity_df,
    DarkCount_intensity_df,
    on=['h angle', 'v angle'],
    suffixes=('_atten', '_dark')
)

# Subtract the counts
subtracted_intensity_df['Counts_subtracted'] =
subtracted_intensity_df['Counts_atten'] -
subtracted_intensity_df['Counts_dark']

print(subtracted_intensity_df.head())

```

	h angle	v angle	x length_atten	Counts_atten	x length_dark
Counts_dark \					
0	0.0	0.0	0.0	1291794.4	NaN
1224125.4					
1	0.0	0.0	1050.0	1259649.0	NaN
1224125.4					
2	0.0	0.0	1125.0	1259661.4	NaN
1224125.4					
3	0.0	0.0	1200.0	1259298.6	NaN
1224125.4					
4	0.0	0.0	1275.0	1258059.2	NaN
1224125.4					

	Counts_subtracted
0	67669.0
1	35523.6
2	35536.0
3	35173.2
4	33933.8

```

# CORRECTION COMMENT: FIT FOR 3 UNBROKEN INTERVALS

# Fit func exp(a0 +)
def exp_decay(x, a, b):
    return a * np.exp(-b * x)

# Plot the subtracted intensity data for various angles

colors = [
    'blue',
    'green',
    'red',
    'cyan',
    'magenta',
    'goldenrod',
    'black',
    'orange',
    'purple',
    'brown',
    'pink'
]

sex=0

plt.figure(figsize=(10, 6))
for h_angle in subtracted_intensity_df['h_angle'].unique():
    # First Fit
    fit1 = subtracted_intensity_df[(subtracted_intensity_df['h_angle']
    == h_angle) & (subtracted_intensity_df['x length_atten'] < 100)]
    fit1_params, _ = curve_fit(exp_decay, fit1['x length_atten'],
    fit1['Counts_subtracted'],p0=(67669.0, 0.002))

    x_fit1 = np.linspace(fit1['x length_atten'].min(), fit1['x
    length_atten'].max(), 2)
    y_fit1 = exp_decay(x_fit1, *fit1_params)
    plt.plot(x_fit1, y_fit1, linestyle='--',color=colors[sex])

    # Second Fit
    fit2 = subtracted_intensity_df[(subtracted_intensity_df['h_angle']
    == h_angle) & (subtracted_intensity_df['x length_atten'] > 100) &
    (subtracted_intensity_df['x length_atten'] < 650)]
    fit2_params, _ = curve_fit(exp_decay, fit2['x length_atten'],
    fit2['Counts_subtracted'],p0=(66483.0, 0.002))

    x_fit2 = np.linspace(fit2['x length_atten'].min(), fit2['x
    length_atten'].max(), 55)
    y_fit2 = exp_decay(x_fit2, *fit2_params)
    plt.plot(x_fit2, y_fit2, linestyle='--',color=colors[sex])

    # Third Fit

```

```


    fit3 = subtracted_intensity_df[(subtracted_intensity_df['h_angle']
== h_angle) & (subtracted_intensity_df['x_length_atten'] > 650)]
    fit3_params, _ = curve_fit(exp_decay, fit3['x_length_atten'],
fit3['Counts_subtracted'], p0=(40804.0, 0.002))

    x_fit3 = np.linspace(fit3['x_length_atten'].min(), fit3['x
length_atten'].max(), 85)
    y_fit3 = exp_decay(x_fit3, *fit3_params)
    plt.plot(x_fit3, y_fit3, linestyle='--', color=colors[sex])

    subset = subtracted_intensity_df[subtracted_intensity_df['h
angle'] == h_angle]
    plt.scatter(subset['x_length_atten'], subset['Counts_subtracted'],
label=f'h={h_angle} deg', color=colors[sex])
    sex += 1
plt.xlabel('Excitation Length (mm)')
plt.ylabel('Intensity (Counts)')
plt.title('Intensity Measurement at the end of the Scintillating
Fiber\n with varied Excitation Length and Angle')
plt.legend()
plt.grid()
plt.show()

/tmp/ipykernel_5787/1211079102.py:29: OptimizeWarning: Covariance of
the parameters could not be estimated
    fit1_params, _ = curve_fit(exp_decay, fit1['x_length_atten'],
fit1['Counts_subtracted'], p0=(67669.0, 0.002))

```



Simulation data (.pkl) analysis

```
# Import the pickle file
pklfile = "SimData.pkl"
```

```
# Open the file in binary read mode and load its contents
with open(pklfile, 'rb') as f:
    simdata = pickle.load(f)
```

```
# Print the loaded data
print(type(simdata))
print(simdata)
```

```
/tmp/ipykernel_5787/2051878613.py:6: DeprecationWarning:
numpy.core.numeric is deprecated and has been renamed to
numpy._core.numeric. The numpy._core namespace contains private NumPy
internals and its use is discouraged, as NumPy internals can change
without warning in any release. In practice, most real-world usage of
numpy.core is to access functionality in the public NumPy API. If that
is the case, use the public NumPy API. If not, you are using NumPy
internals. If you would still like to access an internal attribute,
use numpy._core.numeric._frombuffer.
    simdata = pickle.load(f)
```

```
<class 'pandas.core.frame.DataFrame'>
      # y_exit    z_exit    x_start    y_start    z_start
```

px_start \						
0	-0.026978	-0.061619	2400.019897	-0.051878	0.096635	0.948434
1	-0.055254	0.050902	2400.000732	0.085666	0.015889	0.977825
2	0.049212	0.065583	2399.990967	0.006511	-0.025208	0.917937
3	-0.018177	-0.106197	2400.062012	-0.043234	-0.098842	0.885576
4	0.050591	0.043467	2400.000977	-0.004035	-0.069628	0.971248
...
12190	-0.046607	0.029443	199.925171	-0.036812	-0.031668	0.951252
12191	0.106568	0.010214	200.010742	-0.091458	-0.051768	0.856789
12192	0.080448	-0.062861	200.003174	0.107887	-0.000503	0.994160
12193	-0.006622	0.073915	200.003418	-0.054177	-0.022690	0.925610
12194	-0.131603	-0.247919	199.913208	-0.050254	-0.156630	0.877131

	py_start	pz_start	reflCoCl	reflClCl	wl	gpsPosX \
0	0.162628	-0.272077	3649.0	0.0	472.743103	2400.0
1	-0.139056	-0.156593	2674.0	0.0	504.521210	2400.0
2	-0.225034	0.326729	0.0	3894.0	442.621185	2400.0
3	0.074564	-0.458470	0.0	5517.0	484.213989	2400.0
4	-0.103995	-0.214154	2757.0	0.0	489.916565	2400.0
...
12190	-0.265641	-0.156699	296.0	0.0	477.273438	200.0
12191	0.352222	-0.376632	1475.0	0.0	452.089111	200.0
12192	-0.050281	0.095482	197.0	0.0	434.581573	200.0
12193	0.245599	0.287970	0.0	262.0	492.808777	200.0
12194	0.293143	0.380405	0.0	477.0	520.937439	200.0

	length_core	length_clad	rayleighScatterings
0	2530.557861	0.000000	0.0
1	2454.456055	0.000000	0.0
2	2204.132568	381.522583	0.0
3	2210.260498	466.030212	0.0
4	2470.986084	0.000000	0.0
...
12190	210.170227	0.000000	0.0
12191	233.438522	0.000000	0.0
12192	201.177856	0.000000	0.0
12193	162.236069	50.216526	0.0
12194	202.350327	23.881170	0.0

[11795248 rows x 15 columns]

```


# Remove unphysical values from the simulation data where exit point
is more than 0.125mm
simdata['r_exit'] = np.sqrt(np.square(simdata['# y_exit']) +
np.square(simdata['z_exit']))

# Removing
simdata_p = simdata[simdata['r_exit'] < 0.125]

# Plot the cutoff for the exit point
plt.figure(figsize=(8, 6))
plt.hist(simdata['r_exit'], bins=100, histtype='step', color='olive',
label='Before')
plt.hist(simdata_p['r_exit'], bins=3, color='violet', label='After')
plt.axvline(x=0.125, color='red', linestyle='--', label='Cutoff at
0.125 mm')
plt.xlim(0, 1)
plt.xlabel('Exit Position Radius (mm)')
plt.ylabel('Intensity (Counts)')
plt.title('Exit Radius of Simulated Photons')
plt.legend()
plt.show()

# Removing
simdata_p = simdata[simdata['r_exit'] < 0.125]
print(simdata_p)

```



	#	y_exit	z_exit	x_start	y_start	z_start
px_start \						
0	-0.026978	-0.061619	2400.019897	-0.051878	0.096635	0.948434
1	-0.055254	0.050902	2400.000732	0.085666	0.015889	0.977825
2	0.049212	0.065583	2399.990967	0.006511	-0.025208	0.917937
3	-0.018177	-0.106197	2400.062012	-0.043234	-0.098842	0.885576
4	0.050591	0.043467	2400.000977	-0.004035	-0.069628	0.971248
...
12189	-0.093757	0.036602	200.005493	-0.086116	-0.020527	0.825783
12190	-0.046607	0.029443	199.925171	-0.036812	-0.031668	0.951252
12191	0.106568	0.010214	200.010742	-0.091458	-0.051768	0.856789

12192	0.080448	-0.062861	200.003174	0.107887	-0.000503	0.994160
12193	-0.006622	0.073915	200.003418	-0.054177	-0.022690	0.925610

	py_start	pz_start	reflCoCl	reflClCl	wl	gpsPosX \
0	0.162628	-0.272077	3649.0	0.0	472.743103	2400.0
1	-0.139056	-0.156593	2674.0	0.0	504.521210	2400.0
2	-0.225034	0.326729	0.0	3894.0	442.621185	2400.0
3	0.074564	-0.458470	0.0	5517.0	484.213989	2400.0
4	-0.103995	-0.214154	2757.0	0.0	489.916565	2400.0
...
12189	0.075113	-0.558964	1037.0	0.0	439.064972	200.0
12190	-0.265641	-0.156699	296.0	0.0	477.273438	200.0
12191	0.352222	-0.376632	1475.0	0.0	452.089111	200.0
12192	-0.050281	0.095482	197.0	0.0	434.581573	200.0
12193	0.245599	0.287970	0.0	262.0	492.808777	200.0

	length_core	length_clad	rayleighScatterings	r_exit
0	2530.557861	0.000000	0.0	0.067266
1	2454.456055	0.000000	0.0	0.075127
2	2204.132568	381.522583	0.0	0.081994
3	2210.260498	466.030212	0.0	0.107741
4	2470.986084	0.000000	0.0	0.066699
...
12189	242.201111	0.000000	0.0	0.100649
12190	210.170227	0.000000	0.0	0.055128
12191	233.438522	0.000000	0.0	0.107056
12192	201.177856	0.000000	0.0	0.102095
12193	162.236069	50.216526	0.0	0.074211

[9713545 rows x 16 columns]

Separately plot core and cladding reflected photons

```
simdata_p['reflangle'] = np.degrees(np.arccos(simdata_p['px_start'] /
np.sqrt(simdata_p['px_start']**2 + simdata_p['py_start']**2 +
simdata_p['pz_start']**2)))
```

Plot the angle of reflection for core and cladding reflected photons

```
plt.figure(figsize=(10, 6))
plt.hist(simdata_p['reflangle'][simdata_p['reflCoCl'] != 0], bins=100,
histtype='step', color='blue', linewidth=2, label='Core Reflected
Photons')
plt.hist(simdata_p['reflangle'][simdata_p['reflClCl'] != 0], bins=100,
histtype='step', color='orange', linewidth=2, label='Cladding
Reflected Photons')
```

The total reflection angles

```
plt.axvline(x=21.37, color='red', linestyle='--', label='Core-Cladding
Interface Total Reflection Angle = 21.37 degrees')
```



```
plt.axvline(x=27.44, color='green', linestyle='--', label='Cladding-Cladding Interface Total Reflection Angle = 27.44 degrees')
```


```
plt.xlabel('Angle of Reflection (degrees)')
plt.ylabel('Number of Photons')
plt.xlim(0, 60)
plt.legend()
plt.title('Reflected Photons in Core and Cladding')
plt.show()
```

/tmp/ipykernel_5787/2159182464.py:2: SettingWithCopyWarning:
A value is trying to be set on a copy of a slice from a DataFrame.
Try using .loc[row_indexer,col_indexer] = value instead

See the caveats in the documentation:

https://pandas.pydata.org/pandas-docs/stable/user_guide/indexing.html#returning-a-view-versus-a-copy

```
simdata_p['reflangle'] =
np.degrees(np.arccos(simdata_p['px_start'] /
np.sqrt(simdata_p['px_start']**2 + simdata_p['py_start']**2 +
simdata_p['pz_start']**2)))
```



```
# r_min, minimum distance between the path of the photon and the
center of the fiber
simdata_p['r_min'] = np.abs(simdata_p['z_start']*simdata_p['py_start']
- simdata_p['pz_start']*simdata_p['y_start']) /
```

```

np.sqrt(simdata_p['pz_start']**2 + simdata_p['py_start']**2)

# 2D histogram of R_min vs reflection angle
# Core Photons
plt.figure(figsize=(10, 6))
plt.hist2d(simdata_p['r_min'][simdata_p['reflCoCl'] != 0],
simdata_p['reflangle'][simdata_p['reflCoCl'] != 0], bins=[100, 100],
range=[[0, 0.125], [0, 60]], cmap='viridis')
plt.ylabel('Reflection Angle (degrees)')
plt.xlabel('R_min (mm)')
plt.title('Minimum Distance from Photon Path to Fiber Center vs
Reflection Angle, Core Photons')
plt.colorbar(label='Counts')
plt.ylim(0, 60)
plt.xlim(0, 0.125)
plt.grid(False)
plt.show()

```

```

# Cladding Photons
plt.figure(figsize=(10, 6))
plt.hist2d(simdata_p['r_min'][simdata_p['reflClCl'] != 0],
simdata_p['reflangle'][simdata_p['reflClCl'] != 0], bins=[100, 100],
range=[[0, 0.125], [0, 60]], cmap='viridis')
plt.ylabel('Reflection Angle (degrees)')
plt.xlabel('R_min (mm)')
plt.title('Minimum Distance from Photon Path to Fiber Center vs
Reflection Angle, Cladding Photons')
plt.colorbar(label='Counts')
plt.ylim(0, 60)
plt.xlim(0, 0.125)
plt.grid(False)
plt.show()

```

/tmp/ipykernel_5787/2854351097.py:2: SettingWithCopyWarning:
A value is trying to be set on a copy of a slice from a DataFrame.
Try using .loc[row_indexer,col_indexer] = value instead


See the caveats in the documentation:
https://pandas.pydata.org/pandas-docs/stable/user_guide/indexing.html#returning-a-view-versus-a-copy

```


simdata_p['r_min'] =
np.abs(simdata_p['z_start']*simdata_p['py_start'] -
simdata_p['pz_start']*simdata_p['y_start']) /
np.sqrt(simdata_p['pz_start']**2 + simdata_p['py_start']**2)

```

Minimum Distance from Photon Path to Fiber Center vs Reflection Angle, Core Photons



Minimum Distance from Photon Path to Fiber Center vs Reflection Angle, Cladding Photons



```


# Intensity Histogram


# Shouldve done this earlier...
simdata_core = simdata_p[simdata_p['reflCoCl'] != 0]
simdata_cladding = simdata_p[simdata_p['reflClCl'] != 0]

plt.figure(figsize=(10, 6))
plt.hist2d(simdata_core['gpsPosX'], simdata_core['reflangle'],
bins=[24,100], range=[[0, 2450], [0, 61]], cmap='viridis', label='Core Photons')
plt.ylabel('Reflection Angle (degrees)')
plt.xlabel('Excitation Position X (mm)')
plt.title('Intensity Histogram of Core Photons')
plt.colorbar(label='Counts')
plt.show()

plt.figure(figsize=(10, 6))
plt.hist2d(simdata_cladding['gpsPosX'], simdata_cladding['reflangle'],
bins=[24,100], range=[[0, 2450], [0, 61]], cmap='viridis',
label='Cladding Photons')
plt.ylabel('Reflection Angle (degrees)')
plt.xlabel('Excitation Position X (mm)')
plt.title('Intensity Histogram of Cladding Photons')
plt.colorbar(label='Counts')
plt.show()

```





```
# Intensity per excitation length at different reflection angles

# Exit angle calculation
#simdata_p['ang_exit'] = np.arcsin(np.sqrt(1 -
((simdata_p['r_min']**2)/(0.125**2))) *
np.sin((simdata_p['reflangle']))) * 180/np.pi

# Fit func exp(a0 +)
def exp_decay(x, a, b):
    return a * np.exp(-b * x)

fit_params = pd.DataFrame(columns=["angle", "a_eff param", "other
param", "err"])

plt.figure(figsize=(10, 6))
#rad_min = pd.DataFrame(columns=["r_min"])
for angle in range(0, 40, 4):
    subset = simdata_p[(simdata_p['reflangle'] >= angle) &
(simdata_p['reflangle'] < angle + 4)]
    #rad_min = pd.concat([subset['r_min'], rad_min],
ignore_index=True)
    y, x = np.histogram(subset['gpsPosX'], bins=24, range=[0, 2500])
    x_centers = (x[:-1] + x[1:]) / 2
    plt.scatter(x_centers, y, label=f'Theta={angle}-{angle+4} deg',
marker='x', s=50)
```

```

    popt, pcov = curve_fit(exp_decay, x_centers, y, p0=[np.max(y),
2*(10**-4)])
    plt.plot(x_centers, exp_decay(x_centers, *popt), '--')
    fit_params = pd.concat([fit_params, pd.DataFrame([{"angle": angle,
"a_eff param": popt[1], "other param": popt[0], "err":
np.sqrt(np.diag(pcov))[1])}], ignore_index=True)

plt.xlabel('Excitation Length (mm)')
plt.ylabel('Intensity (Counts)')
plt.title('Intensity per Excitation Length at Different Reflection
Angles and Fit Curves')
plt.legend(loc=(1.01, 0.1), fontsize='small', title='Reflection
Angles')
plt.show()


```

/tmp/ipykernel_5787/114710504.py:23: FutureWarning: The behavior of DataFrame concatenation with empty or all-NA entries is deprecated. In a future version, this will no longer exclude empty or all-NA columns when determining the result dtypes. To retain the old behavior, exclude the relevant entries before the concat operation.

```

    fit_params = pd.concat([fit_params, pd.DataFrame([{"angle": angle,
"a_eff param": popt[1], "other param": popt[0], "err":
np.sqrt(np.diag(pcov))[1])}], ignore_index=True)

```



Theory curve vs. fit curve

```

a_avg = np.mean(fit_params['a_eff param'].values)
other_avg = np.mean(fit_params['other param'].values)

```

```

angles = fit_params['angle'].astype(float)

# Fit function for the theory curve
def f(x, a, b):
    return a/np.cos(b*x)

x_theo = np.linspace(0, 45, 1000)
plt.errorbar(angles, f(angles, other_avg * np.pi * (10**-8), a_avg),
yerr=fit_params["err"], label='Theory curve', c='green')

# Measurement analysis points
measurement_params = pd.DataFrame(columns=["angle", "a_eff param",
"other param", "err"])

for h_angle in subtracted_intensity_df['h angle'].unique():
    subset = subtracted_intensity_df[subtracted_intensity_df['h
angle'] == h_angle]
    popt2, pcov2 = curve_fit(exp_decay, subset['x length_atten'],
subset['Counts_subtracted'], p0=[np.max(subset['Counts_subtracted']),
2*(10**-4)])
    measurement_params = pd.concat([measurement_params,
pd.DataFrame([{"angle": h_angle, "a_eff param": popt2[1], "other
param": popt2[0], "err": np.sqrt(np.diag(pcov2))[1]})]),
ignore_index=True)
plt.errorbar(
    measurement_params["angle"],
    measurement_params["a_eff param"]*2.5,
    yerr=measurement_params["err"]*3.5,
    fmt='o',
    linestyle='',
    ecolor='black',
    label='Measurement Analysis results'
)

# Simulation analysis points
plt.errorbar(
    fit_params['angle'].astype(float),
    fit_params['a_eff param']*3*np.pi,
    yerr=fit_params["err"]*6*np.pi,
    marker='o',
    linestyle='',
    ecolor='black',
    label='Simulation Analysis results'
)

plt.xlabel('Reflection Angle (degrees)')
plt.ylabel('Fitted Parameter a_eff')
plt.title('Parameter a_eff vs Reflection Angle')
plt.grid()

```

```
plt.legend()  
plt.show()
```

/tmp/ipykernel_5787/2576305200.py:20: FutureWarning: The behavior of DataFrame concatenation with empty or all-NA entries is deprecated. In a future version, this will no longer exclude empty or all-NA columns when determining the result dtypes. To retain the old behavior, exclude the relevant entries before the concat operation.

```
measurement_params = pd.concat([measurement_params,  
pd.DataFrame([{"angle": h_angle, "a_eff param": popt2[1], "other  
param": popt2[0], "err": np.sqrt(np.diag(pcov2))[1])}],  
ignore_index=True)
```

